



NeuroAcoustic Transduction Technology

A Personalized, Drug-Free Approach for Insomnia and Sleep Optimization

DELTA BRAINWAVES



ONE INDIVIDUAL, 5 NIGHTS



FIVE INDIVIDUALS, 1 NIGHT

Power of Personalized Sound Sequences (PSS)

Individuals produce a unique pattern of delta brainwaves during the deepest phase of sleep (N3).

Auditory interventions have been shown to:


- Enhance heart rate variability¹
- Reduce cortical hyperarousal²
- Synchronize slow brainwave dynamics³

HOW IT WORKS: *3 Steps to Personalized Sound Therapy*



Step 1: Sleep Data Collection

- Patient wears the Pebble device at home for one night
- Collects unfiltered, 2-channel EEG data (EDF format)



Step 2: Delta Wave Processing

- Detection of delta brainwave activity (0.5–4 Hz) during N3 sleep
- Generation of sound envelopes using Gaussian Fitting
- Creation of 3 personalized MP3 sound sequences





Brainwave Sound Sample

Step 3: The Patient Downloads Their Personalized Sound Sequences (PSS) for Nightly Listening

- Used 30 minutes before bedtime
- Includes a variety of personalized background sounds options tailored to individual preferences

Early Validation in Healthy Volunteers (n = 36):⁴

- 82% reported a relaxation effect
- 67% reported sleep induction
- Randomized surrogate sequences = 0% sleep induction



Insomnia Study (n = 13):⁵


Randomized, single-blind crossover design compared PPS vs. placebo sounds over 3–5 nights

Outcomes:

- Total sleep time: +5%
- REM sleep: +10%
- REM latency: -18%
- Time to fall asleep: up to 40% Reduction

EDF: European Data Format;
EEG: Electroencephalography

References 1. Huwiler S, et al. Effects of auditory sleep modulation approaches on brain oscillatory and cardiovascular dynamics. *Sleep* 2022;45(9):zsac155. 2. Bressler S, et al. A randomized controlled trial of alpha phase-locked auditory stimulation to treat symptoms of sleep onset insomnia. *Sci Rep* 2024;14(1):13039. doi: 10.1038/s41598-024-63385-1. 3. Anderson A. Acoustic enhancement of slow wave sleep—timing and tone: “the details create the big picture”. *Sleep* 2022;10:zsac191 <https://doi.org/10.1093/sleep/zsac191>. 4. Data at myWaves – Computational Neuroscience. 5. Aloulou A, et al. *Sleep Medicine*, in press, 2025 Listening to the sound of your own brain waves enhances sleep quality and quantity.



Join Our Expanding Network of Clinical Researchers

Ongoing studies are investigating NeuroAcoustic Technology for: Chronic insomnia | COMISA (Co-Morbid Insomnia and Sleep apnea) | Sports and military performance | Infant sleep facilitation | Aging-related sleep disturbance